Coming soon! **Next Generation Kendall™ Silicone Bordered 5-Layer Foam Dressing** Clinical effectiveness you demand, quality your patients deserve. Highly breathable film Foam layer Super **Gentle silicone Transfer** absorbent layer

layer

Proven to protect

- Demonstrated effectiveness in the prevention of pressure injuries¹
- Proven to be as absorbent as the industry leaders²
- Seven days of wear-time with our enhanced gentle adhesion
- Thinner, super-absorbent core provides greater absorbency while remaining pliable

Designed for comfort

• Reduces friction and prevents movement in the wound bed

adhesive

- Recommended for use in **neonatal/pediatric and geriatric** populations
- Improved gentle silicone adhesive minimizes trauma to the wound
- Minimizes discomfort to patients by reducing dressing changes
- Unscented, latex-free, and hypoallergenic
- · Variety of sizes to meet the needs of you and your patients

Experience the next generation of wound care.

Get in touch

- 1. Padula WV. Effectiveness and Value of Prophylactic 5-Layer Foam Sacral Dressings to Prevent Hospital-Acquired Pressure Injuries in Acute Care Hospitals: An Observational Cohort Study. J Wound Ostomy Continence Nurs. 2017;44(5):413-419.
- Data on file. Cardinal Health, 2021.



© 2021 Cardinal Health Canada. All Rights Reserved. CARDINAL HEALTH, the Cardinal Health LOGO, ESSENTIAL TO CARE and KENDALL are trademarks of Cardinal Health and may be registered in the US and/or in other countries. All other trademarks are the property of their respective owners. Lit. No. 21-149 -A (11/2021)