

Increased blood flow
generates a
natural healing response

geko™
wound therapy

Wound Therapy

Providing increased blood circulation to promote
wound healing naturally from the inside



The geko™ Wound Therapy device **stimulates** the common peroneal nerve **activating** the calf and foot muscle pumps, increasing venous, arterial and microcirculatory blood flow.^{1, 2, 3}

The blood flow increase prevents and reduces edema⁴ and maintains TCpO₂ – promoting conditions favorable for wound healing.^{3, 5}

For further information on the geko™ Wound Therapy device and the clinical evidence, please visit our website: www.gekowound.ca

Demonstrations / Ordering:

T: +1 888 244 5579

Email: orders@perfusemedtec.com

1. Tucker AT et al. Int J Angiol. 2010 Spring;19 (1): e31-e37
2. Williams KJ et al. Phlebology. 2015 Jun; 30 (5): 365-72
3. Jawad H et al. 2014 Journal Vasc Surg. Vol 2: 160-65
4. Williams KJ et al. Poster. Vascular Society Annual Scientific Meeting, Glasgow November 2014
5. Warwick D et al. Int J Angiol. 2015 April; 34 (2): 158-65

Perfuse
Medtec