



Cutimed® Sorbion®

Handle any level of wound exudate with Sorbion Technology

The smart way to maintain moisture balance

Cutimed Sorbion is designed with Sorbion Technology – its superabsorbent polymers help to maintain moisture balance for any level of exudate², light or heavy, avoiding dehydration³ and maceration⁴. Smart technology, packed into the tiniest details. And it's those tiny details that make a big impact for patients, carers and budget holders.

Discover the benefits, visit www.cutimedcanada.com



¹ Essity group – commissioned laboratory testing at SMTL: Evaluation of absorbency and fluid retention tested against 20 biggest competitors (based on EN 13726). 2024.

² Essity group commissioned laboratory testing at SMTL: Evaluation of absorbency and fluid retention of different superabsorbent wound dressings (based on EN 13726). 2024

³ Cutting K. Made Easy - Cutimed Sorbion Sachet. Wounds International. 2016.

⁴ Cutting KF. Managing wound exudate using a super-absorbent polymer dressing: a 53-patient clinical evaluation. J Wound Care 2009;18(5):200-205.