



Venous Leg Ulcer

Integrated Therapy
Solutions for Best Results

Clean & Debride

Cutimed® Debriclean
Fast and effective wound bed preparation.



Manage Infection & Exudate

Cutimed® Sorbact®
Cutimed® Sorbion®
Hydrofera Blue®
Choose from a wide range of antimicrobial and super absorbent dressings.



Compression Therapy

Compression is the key component for the treatment of venous leg ulceration.



Venous leg ulcer treatment combining moist wound healing and compression therapy works!

A European case study demonstrates the significant benefits of integrating Cutimed® wound care with JOBST® compression therapy for treating venous leg ulcers:

- 85% of cases showed improvement within 12 weeks*
- 53% of cases were completely healed*

*Published in the EWMA Journal, Volume 13, No. 2, October 2013, pp. 19-23.



essity



Cutimed®

HydroferaBLUE®

JOBST®